



GerMedUSA Inc.
THE COMPLETE SHOP FOR SURGICAL INSTRUMENTS

LEEP INSTRUMENTS



WHY SHOULD YOU BUY FROM US?

because


**WE NOT ONLY SELL !
BUT ALSO MANUFACTURE !
OUR OWN SURGICAL INSTRUMENTS**



GerMedUSA, Inc. is one of the few companies who not only sell but manufacture Surgical Instruments. We supply thousands of different types of Surgical Instruments. Our craftsmen pay great attention to detail when producing Medical Surgical Instruments and therefore resulting in quality products.

GerMedUSA supplies a vast number of instrumentation selections for all surgical specialties such as: Cardiovascular, Gynecological, Micro-Surgical, Neurosurgical, Non Scalpel Vasectomy, Plastic Surgery, E.N.T., Dermal, Orthopedic, MIS and Left Handed Surgical Instruments. We also supply Rainbow color Titanium coated Surgical Instruments. All of our Surgical, Veterinary and Dental Instruments have a Lifetime Warranty on them. We also encourage you to send us your surgical instrument ideas or design to possibly help you with manufacturing.

To find out more, please feel free to contact our Customer Service Department or visit us at our website: www.GerMedUSA.com.



Warranty:

All GermedUSA surgical instruments are guaranteed for life against manufacturing defects, provided that the instrument is used for its intended surgical purpose. GermedUSA will either repair or replace the instrument without charge at our discretion. Our liability under this guarantee shall be limited to repair or replacement of defective merchandise.


Our instruments are carefully examined by our surgical instrument experts.

No instrument is shipped unless it meets our standards.

Tungsten carbide instruments are included except for jaw inserts and Gold handles by a separate warranty.

2417 Jericho Turnpike, #333, Garden City Park, NY 11040 (USA)

Tel: 800-330-1322 - 516-358-2180 Fax: 800-260-7177 - 516-358-2182



LATERAL VAGINAL RETRACTOR

G33-220 8-1/4" (21cm), 22 mm x 57 mm blade



TENACULUM FORCEPS

G33-240 BRAUN Tenaculum Forceps 9-3/4"
(24.8 cm), Straight, coated tips



REDUNDANT WALL VAGINAL RETRACTOR

G33-200 7-3/4" (19.5 cm), Medium 22 mm x 57 mm blade

G33-205 7-3/4" (19.5 cm), Large 22 mm x 75 mm blade

GRAVES VIEW-MORE SPECULUM WITH SMOKE EVACUATION TUBE



- G33-170 Medium 4" x 1-1/4" (10.2 cm x 3.2 cm)
- G33-175 Large 4-1/2" x 1-3/8" (11.4 cm x 3.5 cm)

MAXI-VIEW SPECULUM WITH SMOKE EVACUATION TUBE




- G33-180 Pederson Medium 4" x 7/8" (10.2 cm x 2.2 cm)
- G33-185 Pederson Large 4-1/2" x 1" (11.4 cm x 2.5 cm)
- G33-190 Graves Medium 4" x 1-1/4" (10.2 cm x 3.2 cm)
- G33-195 Graves Large 4-1/2" 1-3/8" (11.4 cm x 3.5 cm)

KOGAN ENDOCERVICAL SPECULA



- G33-600 Kogan, 10-3/4" (27.3 cm), narrow tip (3.2 cm x 25 mm), grip-lock ratchet and gauge
- G33-605 Kogan, 10-3/4" (27.3 cm), regular tip (7.2 mm x 25mm), grip-lock ratchet and gauge
- G33-610 Kogan, 10-3/4" (27.3 cm) narrow tip (3.2 mm x 25mm)
- G33-620 Kogan, 10-3/4" (27.3 cm) regular tip (7.2 mm x 25mm)





GRAVES SPECULUM WITH SMOKE EVACUATION TUBE

- G33-120 Small 3" x 7/8" (7.6 cm x 2.2 cm)
- G33-125 Medium 4" x 1-1/4" (10.2 cm x 3.2 cm)
- G33-130 Large 4-1/2" x 1-3/8" (11.4 cm x 3.5 cm)
- G33-135 X-Large 6-3/4" x 1-3/8" (17.1 cm x 3.5 cm)



PEDERSON SPECULUM WITH SMOKE EVACUATION TUBE

- G33-100 Medium 4" x 7/8" (10.2 cm x 2.2 cm)
- G33-105 Large 4-3/4" x 1" (12.1 cm x 2.5 cm)
- G33-110 X-Large 6-3/4" x 1" (17.1 cm x 2.5 cm)



OPEN SIDED SPECULUM (RIGHT) WITH SMOKE EVACUATION TUBE

- G33-15 Pederson Med. 4" x 7/8" (10.2 cm x 2.2 cm)
- G33-155 Pederson Lg. 4-1/2" x 1" (11.4 cm x 2.5 cm)
- G33-160 Graves, Large 4-1/2" x 1-3/8" (11.4 cm x 3.5 cm)
- G33-165 Graves, X-Large 6-3/4" x 1-3/8" (17.1 cm x 3.5 cm)




IRIS HOOKS



G33-400 10" (25.4 cm) long, small, delicate 4 mm hook, with exposed tip.



EMMETT TENACULUM HOOKS



G33-405 Style 1, half curved, small
 G33-410 Style 2, acute angle
 G33-415 Style 3, right angle
 G33-420 Style 4, half curved, large



TISSUE FORCEPS DESCRIPTION



G33-310 8" (20.3 cm) long, 1x2 teeth with exposed tips
 G33-315 8" (20.3 cm) long, 1x2 teeth with coated tips
 G33-320 10" (25.4 cm) long, 1x2 teeth with exposed tips
 31-51 10" (25.4 cm) long, 1x2 teeth with coated tips




DRESSING & TISSUE FORCEPS

G33-300 8" (20.3 cm) long, serrated with coated tips
 G33-305 10" (25.4 cm) long, serrated with coated tips



STERILE DISPOSABLE ELECTRODES

G33-700	Ball (3 mm Diameter)	G33-735	Loop (15 mm width x 10 mm depth)
G33-705	Ball (5 mm Diameter)	G33-740	Loop (20 mm width x 8 mm depth)
G33-710	Square (10 mm width x 8 mm depth)	G33-745	Loop (20 mm width x 10 mm depth)
G33-715	Loop (10 mm width x 8 mm depth)	G33-750	Loop (20 mm width x 12 mm depth)
G33-720	Loop (10 mm width x 10 mm depth)	G33-755	Loop (20 mm width x 15 mm depth)
G33-725	Loop (15 mm width x 5 mm depth)	G33-760	Loop (20 mm width x 20 mm depth)
G33-730	Loop (15 mm width x 8 mm depth)		



GermedUSA Inc.
 THE COMPLETE SHOP FOR SURGICAL INSTRUMENTS

TUNGSTEN CARBIDE Instruments

Why Should you buy from us? Because We not only Sell, but also Manufacture our Surgical Instruments

1617 Garden Turnpike
 #333, Garden City Park
 NY 11040-2814
 Tel: 800-330-1322 / 516-358-2180
 Fax: 800-260-7177 / 516-358-2182
 www.GerMedUSA.com

GermedUSA Inc.
 THE COMPLETE SHOP FOR SURGICAL INSTRUMENTS

20% Off
 SPECIAL OFFER

GermedUSA

LEFT HAND SURGICAL INSTRUMENTS

Why Should you buy from us? Because We not only sell but also Manufacture our own Surgical Instruments

A Complete Shop for Surgical Instruments

GermedUSA Inc.
 THE COMPLETE SHOP FOR SURGICAL INSTRUMENTS

Why Should you buy From us?

WE NOT ONLY SELL BUT ALSO MANUFACTURE our own SURGICAL INSTRUMENTS

SPECIALS

GermedUSA Inc.
 THE COMPLETE SHOP FOR SURGICAL INSTRUMENTS

GENERAL SURGICAL INSTRUMENTS

SCISSORS - FORCEPS - RETRACTORS - NEEDLE HOLDERS

Why should you buy from us? We not only SELL, but also MANUFACTURE our own SURGICAL INSTRUMENTS

GerMedUSA Inc.
 THE COMPLETE SHOP FOR SURGICAL INSTRUMENTS

Tungsten Carbide
SUPER Sharp

@
 www.GerMedUSA.com

Why should you buy from us? we not only Sell but also Manufacture



GermedUSA Inc.
 THE COMPLETE SHOP FOR SURGICAL INSTRUMENTS

SURGERY SETS

- GENERAL SURGERY SETS
- CRANIAL SURGERY SETS
- GENITO - URINARY SURGERY SETS
- EAR, NOSE AND THROAT SURGERY SETS
- OPHTHALMIC SURGERY SETS
- RECTAL SURGERY SETS
- NEUROLOGICAL SET
- ORTHOPEDIC SURGERY SETS
- THORACIC AND CARDIAC SURGERY SETS
- PEDIATRIC SETS

GermedUSA Inc.
 THE COMPLETE SHOP FOR SURGICAL INSTRUMENTS

RAINBOW COLOR TITANIUM COATED SURGICAL INSTRUMENTS

SCISSORS - FORCEPS - RETRACTORS - NEEDLE HOLDERS

Why should you buy from us? Because We not only sell but also manufacture our own surgical Instruments

**2417 Jericho Turnpike, #333,
 Garden City Park, NY 11040 (USA)
 Tel: 800-330-1322 - 516-358-2180
 Fax: 800-260-7177 - 516-358-2182
 E-mail: sales@germedusa.com,
 www.GerMedUSA.com**