

UROLOGY INSTRUMENTS



*GerMed*USA™

Address: 2417 Jericho Turnpike, #333, Garden City Park, NY 11040 (USA).
Phone: +1(516) 997-2180 Fax: +1(516) 358 2182 Email: sales@germedusa.com

About the Company	2
Autoclavable Cystoscopes	3
Resectoscope - Working Elements	4
Resectoscopy Sheaths Obturators, Working Elements	5
Electrodes	6
Cystoscope Sheaths	7
Optical Instruments	8
Flexible Instruments	9
Urethrotome	10

Special models available upon request

ABOUT US

GerMedUSA Inc. is the leading manufacturer of high-quality surgical instruments worldwide. At most economical prices, we offer premium quality surgical instruments, immediately available for delivery anywhere in the world. Our broad range of surgical instruments undergoes a series of quality control testing to get excellent products. We offer a wide selection of Surgical Instruments, including; Dental Instruments, Orthopedic Surgical Instruments, Operating Room Surgical Instruments, Podiatry Surgical Instruments, Rainbow Surgical Instruments, Left Hand Instruments, Thoracoscopy, Gynecological Instruments, Dermal Instruments, Ear Instruments, Eye Instruments, Hemostatic Forceps, Cardiovascular Thoracic Instruments, Floor Grade Economy Instruments, and Tungsten Carbide Insert Needle Holders and Scissors.

GerMedUSA Inc. journey that began three decades ago has successfully become a symbol of Excellence, Innovation, and Reliability for veterinarians worldwide. Strictly following the **FDA** guidelines and being **ISO 13485 certified**, all our procedures are documented, and every process is being monitored at every level. We have acquired an unmatched and unique understanding of manufacturing. Our state-of-the-art manufacturing facilities and years of experience have helped us keep prices low and provide quality products. We believe more in building a strong relationship than just selling the instruments. Therefore, we enjoy providing superior customer services from the first query to the last and offer free sharpening on all our Scissors for the life of the instrument & all Dental Elevators/Luxating Elevators.

Yet, **GerMedUSA Inc.** believes that there is no end when it comes to success and perfection. Therefore, we continuously work with top-notch professionals to understand the changing market needs and future trends to develop the instruments that make them fall in love with their job. Our team of Expert Consultants is available for you whenever you need them. Order a free CD-Rom "Our electronic catalog" to save your time and deliver your order without any hassle.

Want to talk to us? Start a live chat or reach out to us via social media platforms. Or if you are looking for an instrument that we do not have, please get in touch, we will manufacture it for you.


We would love your honest feedback about what we are doing great and what improvements we need. Thanks in advance for your time and input.

Guarantee

All our surgical instruments are guaranteed to be free from any defects in material and workmanship. Any material or instrument found to be defective will either be replaced or repaired based on the evaluation. This service is free of any cost. Our responsibility under guarantee shall be limited to only replacement and repaired policy of defective merchandise.

All telescopes are equipped with a 3-way variable cable connector to connect the following systems:

ACM, BRITISH, WOLF, HSW, STORTZ, OLYMPUS




Cystoscope– Telescopes, length 302 mm, autoclavable at 2.3 bar (134°)

71-00094	4.0 mm, 0°	\$5,090.00
71-00095	4.0 mm, 12°	\$5,230.00
71-00096	4.0 mm, 30°	\$5,230.00
71-00097	4.0 mm, 70°	\$5,470.00




Fiber Light Cables, Autoclavable at 2 bar (134°C)

Endoskop-Adapter		180 length		230 length		300 length	
Rema / Storz	4.8	71-00200	\$225.00	71-00210	\$220.00	71-00220	\$700.00
	6.0	71-00201	\$475.00	71-00211	\$575.00	71-00221	\$1,310.00
Winter / Ibe	4.8	71-00202	\$280.00	71-00212	\$330.00	71-00222	\$790.00
	6.0	71-00203	\$475.00	71-00213	\$585.00	71-00223	\$1,320.00
Wolf / HSW	4.8	71-00204	\$285.00	71-00214	\$345.00	71-00224	\$800.00
	6.0	71-00205	\$480.00	71-00215	\$575.00	71-00225	\$1,330.00
ACM / British	4.8	71-00206	\$275.00	71-00216	\$345.00	71-00226	\$800.00
	6.0	71-00207	\$480.00	71-00217	\$575.00	71-00227	\$1,330.00



74-02100	\$2000.00
----------	-----------

Working element **PASSIVE**, electrode
in resting position inside the sheath




74-02110	\$2000.00
----------	-----------

Working element **ACTIVE**, electrode
in resting position outside the sheath



74-02120	\$1,600.00
----------	------------

**Continuous flow sheath 24 CH fixed
(complete with inner sheath)**



74-02121	\$1,600.00
----------	------------


**Continuous flow sheath 24 CH rotating
(complete with inner sheath)**

74-02114	\$960.00
----------	----------

**Resectoscope Sheath 24 CH
(without picture) with one fixed stopcock**

74-02115	\$1,050.00
----------	------------

**Resectoscope Sheath 24 CH
(without picture) with two fixed stopcocks**



74-02117	\$240.00
----------	----------

**Standard Obturator 24 CH for
Resectoscope Sheath**



74-02130	\$520.00
----------	----------

**Deflecting Obturator 24 CH for
Resectoscope Sheath**



74-02135	\$580.00
----------	----------

**Visual Obturator 24 CH for Resectoscope
Sheath**



74-02030	\$1,200.00
----------	------------

**Working element with Albarran deflector,
1 instrument channel**

74-02031	\$1,260.00
----------	------------

**Working element with Albarran deflector,
2 instrument channels**



74-02200	Loop electrode, curved 24 CH	\$120.00
74-02201	Loop electrode, curved 27 CH	\$120.00



74-02202	Loop electrode, straight 24 CH	\$120.00
----------	--------------------------------	----------



74-02203	Knife electrode 24 CH	\$120.00
----------	-----------------------	----------



74-02204	Ball electrode 24 CH 3 mm	\$120.00
----------	---------------------------	----------



74-02205	Ball electrode 24 CH 5 mm	\$120.00
----------	---------------------------	----------



74-02206	Roller electrode 24 CH 3 mm	\$120.00
----------	-----------------------------	----------



74-02207	Roller electrode 24 CH 5 mm	\$120.00
----------	-----------------------------	----------



74-02208	Conical electrode 24 CH	\$120.00
----------	-------------------------	----------



74-02217	Vaporization electrode 3 mm	\$120.00
----------	-----------------------------	----------



74-02218	Vaporization electrode 5 mm	\$120.00
----------	-----------------------------	----------



74-02215	Spike electrode 3 mm	Call for code
----------	----------------------	---------------



74-02216	Spike electrode 5 mm	Call for code
----------	----------------------	---------------

High-Frequency Cable for use with Resectoscopy - Working element, 3m long.




74-02250	HF-cable for Aesculap, Berchtold, Martin, 4 mm	\$60.00
----------	--	---------



74-02251	HF-cable for Storz, Erbe, 5 mm	\$60.00
----------	--------------------------------	---------




74-02252	HF-cable for Valleylab, Bowie, 8 mm	\$60.00
----------	-------------------------------------	---------




Cystoscope sheaths with 2 fixed stopcocks including standard obturator


74-02010	17 CH	\$1,000.00
74-02011	19 CH	\$1,000.00
74-02012	21 CH	\$1,000.00
74-02013	23 CH	\$1,000.00
74-02014	25 CH	\$1,000.00




74-02020	Diagnostic adapter-bridge	\$320.00
----------	---------------------------	----------



74-02021	Adapter-bridge with one channel	\$410.00
----------	---------------------------------	----------



74-02022	Adapter-bridge with two channels	\$480.00
----------	----------------------------------	----------



74-02050	Biopsy forceps, double-action jaws	\$1,070.00
----------	------------------------------------	------------



74-02053	Biopsy forceps, double-action jaws	\$1,070.00
----------	------------------------------------	------------




74-02052	Hook scissors, double-action	\$1,070.00
----------	------------------------------	------------



74-02054	Biopsy punch, double-action jaws	\$1,130.00
----------	----------------------------------	------------



74-02051	Stone grasping forceps, double action jaws	\$1,070.00
----------	--	------------



74-02055	Stone grasping forceps for use with 25CH sheath	\$1,600.00
----------	---	------------






Button Electrodes flexible, 53 cm		
74-02085	5 CH	\$120.00
74-02086	7 CH	\$120.00




Needle Electrodes, flexible, 53 cm		
74-02088	5 CH	\$120.00
74-02089	7 CH	\$120.00




Loop Electrodes, flexible, 53 cm		
74-02090	5 CH	\$120.00
74-02091	7 CH	\$120.00



Scissors, flexible				
Length	40 cm		60 cm	
5 CH	74-00612	\$295.00	74-00613	\$630.00



Foreign body forceps, flexible				
Length	40 cm		60 cm	
5 CH	74-00615	\$295.00	74-00616	\$630.00
7 CH	74-00601	\$295.00	74-00617	\$630.00




Biopsy forceps, flexible				
Length	40 cm		60 cm	
5 CH	74-00618	\$295.00	74-00619	\$630.00
7 CH	74-00610	\$295.00	74-00620	\$630.00




Stone basket, flexible		
Length	60 cm	
5 CH	74-00608	\$760.00

(Also available as half-flexible or with insulated spiral)



74-02100	Working element PASSIVE, Knife in resting position inside the Sheath	\$2,000.00
----------	--	------------




74-02400	Urethrotome Sheath with Instrument Channel	\$920.00
----------	--	----------



74-02401	Obturator for Urethrotome Sheath, 21 CH	\$270.00
----------	---	----------



74-02403	Guide Tube for insertion of a Ballon Catheter	\$490.00
----------	---	----------




74-02402	Urethrotome Bridge, 1 channel, 5 CH	\$420.00
----------	-------------------------------------	----------



74-02410	Urethrotome Knife, Straight, Cold Cutting	\$130.00
----------	---	----------



74-02411	Urethrotome Knife, Round, Cold Cutting	\$130.00
----------	--	----------



74-02500	Evacuator (Ellik)	\$250.00
----------	-------------------	----------



74-02530	Syringe 150 ml	\$350.00
----------	----------------	----------